



**NovaTec® premium 15-3-20 (+3 +10)**

(83017834)

	3,40
	1,10
calibration	1,24
throw wide	16,2

**Attention!! The data given in this setting chart is based on laboratory tests. An additional cross-check on a mobile test kit is recommended.**

ZA-TS / ZG-TS	mm		kg	%	mm	kg	%	mm	kg	%	mm	kg	%	mm	kg	%
	1	2														
TS 1	15	37	600	A	1	500	1	21	500	1	31	500	19	-4	170	
	16	37	600	A	1	500	1	20	500	1	30	500	19	-3	170	
	18	41	600	A	1	720	1	19	600	1	29	500	19	-1	209	
	20	45	600	A	1	720	1	17	600	1	27	500	19	1	249	
	21	45	720	A	1	600	1	18	600	1	28	500	22	-1	215	
TS 2	21	11	600	B	3	720	3	26	600	3	36	500	24	-3	170	
	24	14	600	B	3	720	3	25	600	3	35	500	25	0	193	
	27	19	600	D	X	500	1	23	600	1	33	500	25	1	247	
	28	19	720	B	3	800	3	23	720	3	33	600	28	0	203	
	30	21	800	B	3	900	3	23	800	3	33	720	30	0	204	
TS 3	20	11	600	B	3	720	3	27	600	3	37	500	24	-4	170	
	24	12	600	C	2	720	2	24	600	2	34	500	28	-3	170	
	27	15	600	C	3	720	3	21	600	3	31	500	28	-1	199	
	28	16	600	C	3	720	3	21	600	3	31	500	28	0	214	
	30	19	600	C	3	720	3	20	600	3	30	500	29	0	244	
	32	21	720	C	3	800	3	20	720	3	30	600	32	0	218	
	33	21	800	C	3	900	3	21	800	3	31	720	34	-1	208	
	36	24	900	C	3	900	3	21	800	3	31	720	36	-1	224	
39	27	900	D	X	720	1	19	800	1	29	720	37	0	258		
40	28	900	D	X	800	1	19	900	1	29	800	37	1	269		